PULPDENT

ACTIVA™ Presto™
Universal Stackable Composite

Light Cure • Fluoride Releasing • Moisture Friendly • Radiopaque

XP-VPF-IN-01w
REV: 10/2019

ACTIVA™ PRESTO™

PRODUCT DESCRIPTION

ACTIVA PRESTO™ is an esthetic, light cure material that contains calcium, phosphate, and fluoride in a durable, wear-resistant resin matrix indicated for all restorative procedures. The material is stackable, does not slump, and holds its shape prior to light curing.

ACTIVA PRESTO is tougher than traditional composites. It contains a patented, shock-absorbing, rubberized-resin component that resists chipping and fracture. This increases flexural strength while maintaining the esthetics, hardness, compressive strength, and durability required of dental composites.

ACTIVA PRESTO is a urethane-based resin that contains no Bisphenol A, no Bis-GMA, and no BPA derivatives.

INDICATIONS

Recommended as a filling material for pits, root surface cavities, and Class I, II, III, IV and V restorations where there is no pulpal involvement.

CONTRAINDICATIONS

Not indicated for direct placement on the exposed pulp.

INSTRUCTIONS FOR USE

1. Isolate and prepare tooth to receive a restoration. Bevel and round off enamel margins. In Class V lesions, bevel or undercut enamel.
2. Place appropriate pulp protection, if indicated.
3. Etch and apply a bonding agent following the manufacturer's instructions for use.
4. **ACTIVA PRESTO** is a dentin and enamel replacement material. Place a 19-gauge applicator tip on the syringe and apply in 2mm increments, light curing for 20 seconds between each layer.
5. Finish and polish as usual.

STORAGE AND HANDLING

- Store tightly sealed in original container at cool room temperature. Avoid direct light, extremes of temperature, contamination and sources of ignition.
- Shelf life of unopened product: 2 years from date of manufacture.
- Re-cap immediately after use.

CAUTION

Uncured material may cause eye or skin irritation on contact. Dental professionals should wear safety glasses and surgical gloves. Do not exceed manufacturer's recommended curing time for the light you are using.

Note: Multi-dose syringes should be encased in a fresh protective barrier for each patient. For additional information, refer to: <http://www.fda.gov/Medical Devices/Products and Medical Procedures/Dental Products/ucm4044 72.htm>.

PHYSICAL PROPERTIES

Light cure setting time:	20 ± 2 seconds
Depth of light cure:	2.6mm
Percentage filler by weight:	70%
Fluoride release 1 day:	230 ppm
Fluoride release 7 days (cumulative):	940 ppm
Flexural strength:	108 MPa ± 5 MPa / 15,660 ± 725 psi
Flexural modulus:	8.0 GPa
Compressive strength:	314 MPa ± 8 MPa / 45,530 psi ± 1160 psi
Diametral tensile strength:	53 MPa ± 2 MPa / 7,685 psi ± 290 psi
Specific gravity	1.820 ± 0.02
Deflection at break	0.8mm
Water sorption (1 week):	1.6%
Solubility	0.08%
Polymerization shrinkage:	2.3%
Radiopacity	2.5 AI (250%)

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